

## PRODUCT PROFILE

**DentaLife**  
AUSTRALIA

# Porcelain Etch Gel

<b>Composition</b>	Hydrofluoric Acid 9.6% w/w in stable, water soluble gel.
<b>Uses</b>	Porcelain Etch Gel has the ability to rapidly, chemically etch porcelain and ceramic restorations to facilitate resin bonding procedures.
<b>Advantages</b>	<ul style="list-style-type: none"><li>• Effective, rapid etching in 45 to 60 seconds</li><li>• Thixotropic gel formula for controlled and safe handling</li><li>• Simple syringe dispensing system</li></ul>
<b>Contraindications</b>	Intraoral use in difficult or uncooperative patients.
<b>Precautions</b>	Avoid contact with skin and do not swallow. Procedure should involve the use of neutralising agent. Clean up excess gel or spills immediately.
<b>Clinical Procedure</b>	Isolate the restoration intraorally with rubber dam and/or barrier paste. Apply a small amount of Porcelain Etch Gel over the porcelain or ceramic requiring etching. Leave for 45 to 60 seconds. Evacuate from the mouth with water and high volume suction immediately after etching time is completed. Neutralise gel where possible.
<b>Poison Schedule</b>	Not Applicable
<b>First Aid</b>	For advice, contact a Poisons Information Centre (Australia) 13 11 26, (Thailand) +66 (0)2 547 0100
<b>Regulatory</b>	Australia : Exempt Thailand : Manufacturing License No. Sor Phor 15/2558
<b>Presentations</b>	2.5 mL Luer Lock Syringe in presentation pouch Product Code : DL1511
<b>Storage</b>	Store below 25°C in a cool, dry place
<b>Shelf Life</b>	2 years from date of manufacture as indicated on pack.

Developed in Australia by : Dentalife Australia Pty Ltd

[www.dentalife.com.au](http://www.dentalife.com.au)

Manufactured by : Pharmadesign Co., Ltd.

199/3 Moo 21, Bangpleeyai Bangplee Samutprakarn Thailand 10540

PD Version 1 Issued: 04/2021