

## PRODUCT PROFILE

**DentaLife**  
AUSTRALIA

# HYPOCHLOR

## 4% FORTE Solution

Composition	Sodium Hypochlorite 45g/L Equivalent to 4.5% w/v available Chlorine. At expiry date: Sodium Hypochlorite 20 g/L. Equivalent to 2.0% w/v available Chlorine.
Uses	Purpose designed, irrigation solution containing Sodium Hypochlorite for use in endodontic procedures to assist debridement and cleaning of the root canal. Also aids in chemical breakdown of pulp soft tissue.
Advantages	<ul style="list-style-type: none"><li>• Registered TGA product for endodontic irrigation manufactured under Therapeutic licence in Australia (replaces non-approved solutions).</li><li>• Formulated specifically for use in dental endodontic procedures</li><li>• Contains pharmaceutical wetting agent to support the specific function of dental irrigation activity and response</li><li>• Formulation is standardised and stabilised</li><li>• Refer Aust. Dental Journal 2006;51:(3):245-251, "Dissolution of porcine incisor pulps in sodium hypochlorite solutions of varying compositions and concentrations" for supporting information</li></ul>
Contraindications	Not for subcutaneous injection.
Precautions	Do not irrigate beyond apex of canal. Do not swallow.
Clinical Procedure	Transfer solution to a disposable syringe fitted with a small gauge needle. Carefully insert the needle into prepared root canal until slight resistance is felt then withdraw about 2mm. Express ENDOSURE HYPOCHLOR 4% Forte solution slowly into the canal to assist debridement and cleaning of the root canal and aid chemical breakdown of pulp soft tissue. DO NOT bind the needle or inject beyond the apex or root canal.
Poison Schedule	Not scheduled
ARTG Listing	160242
First Aid	If poisoning occurs, contact a doctor or Poisons Information Centre. Phone 13 11 26.
Presentations	500mL Dispenser Bottle      Product Code : DL0615 1.25L Cartridge Refill Bottle      Product Code : DL0611
Storage	Store below 25 °C and away from direct sunlight.
Shelf Life	18 months from date of manufacture, indicated on bottle. Shelf life is reduced if cap is not returned immediately after dispensing.



Developed and Manufactured in Australia by: **Dentalife Australia Pty Ltd**  
9/505 Maroondah Hwy Ringwood, Vic 3134  
[www.dentalife.com.au](http://www.dentalife.com.au)