



# CHLORHEXIDINE 10% Solution

<b>Composition</b>	Chlorhexidine Gluconate BP 100mg/mL (10% w/v), Neutral pH
<b>Uses</b>	Chlorhexidine solution, no additives, for professional dental use. Also used in dentistry for cleansing and sanitisation of dental impressions.
<b>Advantages</b>	<ul style="list-style-type: none"> <li>• Highly surfactant solution (strong detergency)</li> <li>• Solution contains only Chlorhexidine and no other excipients ie. NO flavours, sweeteners, preservatives, colourants (dyes) or other excipients</li> <li>• Pharmaceutical grade Chlorhexidine in high purity water</li> </ul>
<b>Contraindications</b>	Not for subcutaneous injection.
<b>Precautions</b>	May cause hypersensitivity in some patients. Irritant, avoid contact with the eyes.
<b>Clinical Procedure</b>	<p>Chlorhexidine 10% Solution is a concentrate solution for a broad range of applications in the dental environment. The 10% solution is a concentrate which can be diluted to a concentration for the desired application. Use distilled or purified water for dilution. Do not add other chemicals or detergents as these may react with the chlorhexidine.</p> <p>Check the industry guidelines, procedural requirements or recommended concentration of Chlorhexidine for intended use before preparing a dilution. Use the measuring cup provided to prepare the ratio of 10% Chlorhexidine in Purified Water to dilute to the required concentration. Do not store diluted solutions. Dispose of diluted solutions at the end of use.</p>
<b>Poison Schedule</b>	S7
<b>First Aid</b>	If poisoning occurs, contact a doctor or Poisons Information Centre. Phone 13 11 26.
<b>Presentation</b>	200mL Amber Glass Dispenser Bottle                      Product Code : DL5210
<b>Shelf-life</b>	2 years from date of manufacture, indicated on bottle
<b>Storage</b>	Store below 25°C