

Clinicare

Suction Cleaner

Daily & Weekly



For the daily and weekly cleaning and maintenance of dental aspiration and suction systems.

Liquid concentrate for simple preparation of solution for daily and weekly flushing of

- Aspiration systems
- Suction System
- Spittoons

Advantages:

- Three actions
 - decontaminates, sanitises and deodorises
- Alkaline cleaner for removal of daily and weekly deposits
- Simple and safe application
- Non disruptive foam
- Economical concentrate for preparation of:
 - ~ 2% dilution for daily cleaning and deodorising
 - ~ 6% dilution for weekly maintenance and decontaminates
- Aldehyde free, Phenol free, Chlorine free



Combined with the daily and weekly cleaning of your dental suction, ensure you schedule a monthly recovery (shock) treatment for your treatment centre. This will further assist in the maintenance of your evacuation system by further sanitising and removing debris.

Dentalife
AUSTRALIA

Manufactured by:

Dentalife Australia Pty Ltd
Factory 9 /505 Maroondah Highway
Ringwood, Victoria 3134 Australia
Phone (03) 9879 1226
www.dentalife.com.au

