

Clinicare Suction Cleaner

Monthly Recovery



Combined with the daily and weekly cleaning of your dental suction, ensure you schedule a monthly recovery (shock) treatment for your treatment centre. This will further assist in the maintenance of your evacuation system by further sanitising and removing debris.

Liquid concentrate for simple preparation of solution for monthly flushing of

- **Aspiration systems**
- **Suction System**
- **Spittoons**

Advantages:

- Three action
 - sanitises, decontaminates and deodorises
- Acidic cleaner for removal of build up
- Simple and safe application
- Non disruptive foam solution
- Does not contain: Aldehydes, Phenols or Chlorine

Reorder code:
DL2741 - 1 L Chamber Bottle



For the daily and weekly cleaning and maintenance of dental aspiration and suction systems.

Dentalife
AUSTRALIA

Manufactured by:

Dentalife Australia Pty Ltd
Factory 9 /505 Maroondah Highway
Ringwood, Victoria 3134 Australia
Phone (03) 9879 1226
www.dentalife.com.au

