


# ActivEtch 37% Gel & Solution


The ideal etch for  
improved bond strengths



- Ideal etch for both enamel and dentine surfaces
- Viscous solution that is easy to place and stays where you want, without slumping
- Self-limiting etchant, giving you the peace of mind of not over etching the surface
- Blue in colour to allow for ease of placement
- Easy to rinse off with water leaving no residual
- Australian Made




## Porcelain Etch



- Viscous 9.6% Hydrofluoric acid
- For use on porcelain veneers, inlays, onlays and for porcelain repairs
- Easy to place and control etch
- Effective etching in 45-60 seconds
- Simple syringe dispensing system
- Use with Dentalife Silane

## Silane



- Surface modifying agent
- For use after porcelain etch
- Promotes the bond strength between reins, porcelain and ceramic surfaces



**Manufactured by:**  
Dentalife Australia Pty Ltd  
Factory 9 /505 Maroondah Highway  
Ringwood, Victoria 3134 Australia  
Phone (03) 9879 1226  
www.dentalife.com.au